



# BLOCK WARMER

Suitable for heating sample in the test tube



## FEATURES

- Error code display
- Fixed safety circuit
- LED display for easy reading
- Adjustable temperature is up to 150°C
- High temperature stability is up to + 0.2°C
- Adjustable timing range of 1 min to 99h59min
- 3D digital wheel adjustment of time and temperature
- Automatic calibration with external PT1000 temperature sensor
- Machined aluminum heating unit allows for highly precise control of temperature
- Included PT1000 temperature sensor allows temperature control directly in the sample vessel

 **ADDRESS**  
PULSE SCIENCE CO., LTD.  
28/9 MOO 1 BANGSRIMUANG, MUANG,  
NONTABURI 11000 THAILAND

 **TELEPHONE & FAX**  
+66 2886 7808  
+66 2886 7809

 [FACEBOOK.COM/PULSESCIENCEGROUP](https://www.facebook.com/pulsesciencegroup)

 [LINKEDIN.COM/COMPANY/PULSESCIENCE](https://www.linkedin.com/company/pulsescience)

 [WWW.PULSESCIENCE.CO.TH](http://www.pulsescience.co.th)



URS is a member of Registrar of Standards(Holdings)Ltd.

