

PULSE SCIENCE

OxiWell

Hyperbaric CHAMBER




The *Healing*
Power with Oxygen **100%**



WellSeat Hyperbaric Chamber

Benefits of Hyperbaric Oxygen



Size

L1500 * W600 * H1170 mm

Pressure

1.3 ata

Soft sitting chamber has one adjustable sofa chair inside, sitting in a longer, higher, wider oxygen chamber, or adjust sofa chair to lying down. More than health benefits, we like to bring comfort to users also.



Delay aging

The proven test results show that hyperbaric oxygen can prolong cell telomerase and reduce senescent cells.



Improve memory

Improve blood supply to the brain, strengthen the cerebral cortex, improve the vitality of brain cells, and improve learning efficiency.



Improve physical fitness

Regular micro-hyperbaric oxygen inhalation can improve subhealth, increase cell vitality, and enhance physical function.



Improve hypoxia

It can improve the relative hypoxia caused by drinking, smoking, lack of sleep time, etc



Relieve fatigue

Accelerate the proliferation of fibroblasts and endothelial cells, speed up the recovery of the body, and improve vitality.



Improve sexual function

Adrenal blood supply and sufficient oxygen under hyperbaric oxygen improve the body's hypoxic state and maintain vigorous energy.



Beauty

Strengthen the aerobic metabolism of cells, reduce melanin pigmentation, reduce ecchymosis, and beautify the skin.



Improve vision and presbyopia

Improve retinal ischemia and hypoxia, promote tissue metabolism, and clear retinal vein blockage. Improve retinal ischemia and hypoxia, promote tissue metabolism, and clear retinal vein blockage.



Air Compressor

Capacity	80~120lpm
Power of compressor	480w
Power of cooler	390w
Timer	HBOT time setting



Oxygen Concentrator

Oxygen flow rate	5/10lpm
Power	550w
Oxygen purity	90~96%
Warning	Power fail, Purity, Pressure warning

WellPro

Hyperbaric Chamber



Size

Ø760 * 2160 mm

Pressure

1.5 ata

Thickness

4mm aluminum alloy



Refrigeration Dryer

water drainage & cooling

Capacity	420lpm
Power	390w
Size	26.7*18*20mm
Weight	7kgs
Working pressure	0.8~1.0MPa




Oxygen Concentrator

2 in 1 oxygen concentrator

Oxygen flow rate	5lpm
Power	700w
Oxygen purity	90~96%
Noise level	55dB
Compression pressure	0~200kpa
Compression flow rate	200kpa,75lpm
Touch screen	7inch LCD screen
Oxygen output pressure	0.06MPa


Hard lying chamber 1.5ATA means the housing need to beary higher pressure extrusion, compare to 1.3ata hard oxygen chamber, we thickened the material and add more fasten fittings on it also. Except the inside latex mattress for wonderful experience, Side walk-in big window is more friendly for our users. Bottom wheels make oxygen chamber move easier.



HBOT model option

Operate panel


2 in 1 oxygen concentrator integrated oxygen concentrator and air compressor into 1 package. Its small footprint takes smaller installation room and easier connection. This system with 7 inch lcd touch screen, and built-in oxygen therapy options and HBOT knowledge.





WellSoft

Hyperbaric Chamber



Size

Ø800 * 2250 mm

Pressure

1.5 ata

For higher HBOT pressure, the more durable oxygen chamber material should be update the same. Not only change on pressure valves, oxygen concentrator & cooling compressors, also material of oxygen chamber. Compare to 1.3ata O₂ chamber, 30% more thickness on the 1.5ata oxygen chamber TPU material.



Air Compressor

Capacity	80~120lpm
Power of compressor	480w
Power of cooler	390w
Timer	HBOT time setting



Oxygen Concentrator

Oxygen flow rate	5/10lpm
Power	550w
Oxygen purity	90~96%
Warning	Power fail, Purity, Pressure warning

PULSE SCIENCE

HEAD OFFICE ADDRESS

PULSE SCIENCE CO.,LTD.
28/9 Moo 1 BangSriMuang, Muang,
Nonthaburi, 11000 Thailand


FACTORY ADDRESS

PULSE SCIENCE CO.,LTD.
8/8 Moo 8 LamPho, BangBuaThong,
Nonthaburi, 11110 Thailand

PHONE


+66 2886 7808
+66 2077 7377
+66 86 971 9067



 Pulse Science Group

 Pulse Science

 Pulse Science

 www.pulsescience.co.th

EMAIL SALE@PULSESCIENCE.CO.TH