

# BIOSAFTEY MOBILE UNIT **BMU-TRAILER**



“ Combining medical innovation with state-of-the-art technology is putting the best result on wheels and certified it with the highest qualifications. ”



Sample storage



Pressure Control



Sink and Water system

It's a sample collection cabinet on wheels. It is designed to meet the Centers for Disease Control (CDC) standards for operator safety.

## FEATURES

- Complete separation between medical staff and patient
- Operates under positive pressure for medical staff protection.
- The air will flow from the top of booth and go through HEPA (H 13) filter to provide clean air environment inside.
- The air condition to ensure proper temperature inside the room.
- UV-C or Ozone for sterilization of the internal
- Integrated disinfection spray of sterilization of outside
- Integrated with sample storage system for preserve the sample at 2-8°C