

BIOSAFETY

ISOLATOR CLASS III

Biosafety Isolator Class III is fully enclosed cabinet designed for the protection of personnel or technician and working environment from various microbial contamination while working. Suitable for microbiological applications such as viruses, bacteria and cell or tissue culture.




Feature

- Negative pressure control system
- Inside the cabinet is an air cleaner with HEPA filter, the best grade filter, which can capture 99.99% of small particles as virus.
- The UV-C light disinfection installed makes sure the virus loses its disease-causing ability.
- Fully enclosed cabinet providing a complete barrier or front glass with glove port enclosure and sleeve gloves.



Glove port



Pass box