

WORKSTATIONS

EnviSafe cell

is a Biosafety cabinet Class II, design for visually examining in the clean area. Contaminated area surrounded with negative pressure to fully prevent leakage of contaminated air.




FEATURES

- Semi-closed Biosafety cabinet Class II
- Table top is made of stainless steel integrated with heating areas to give an optimal working condition
- High resolution touch-screen display
- Multiple protection design increases security level significantly
- Innovative surrounded air intake grill design
- ECM brushless DC motor with good energy efficiency
- Supply/Exhaust HEPA filter 99.99% efficiency @ 0.3 micron (Class H14)
- 254 nm UV lamp with 6mm anti-ultraviolet tempered glass
- 1,000 Lux fluorescent lamp
- Ergonomic 10 degree angle

STANDARD BSC CLASS 2



Stem Cell Isolation or Cell Handling is ideally conducted in a clean environment.



Innovative surrounded air intake grille design. All contaminated air in the work area is enclosed inside the cabinet.