

 **ecopipette**
SINGLE CHANNEL PIPETTES

ecopipette Innovation

The **ecopipette** is packed using recyclable, biodegradable materials and constructed from some of the most renewable durable resources available, to ensure an extended lifetime.
Capp ecotrade™ program: Return any brand of old pipettes in exchange for discounted new ecopipettes™.
At Capp ApS, we will ensure an environmentally safe waste discharge of all old pipettes.
Ask your favourite authorized Capp distributor for more information about our ecotrade™ program.

Capp ecopipette™ Unique Advantages:

- Lockable volume controller requiring 50% fewer turns
 - Display is fully visible during hand operations
- Reduced force requirement for aspiration and blow-out
- The only pipette with fixed and variable exchangeable counter mechanism
 - Large pad for ease of tip ejection
- Fully autoclavable pipette body (no disassembly required)
 - Double o-rings for optimal tip seal
- Compatible with all major pipette tip brands
 - Robust design includes metal tip base
- Lowest incidence of breakage resulting from accidental dropping or tips mounting



Ordering Information

ecopipette Starter Packs

Cat. No.	Description
SKE-03-XX	3 pipettes free of choice, stand, 3 racks of tips and a timer



ecopipette variable volume

Cat. No	Description	Inaccuracy %	Imprecision %	Colour
C02-1S	Capp pipette, variable vol. 0.2-2 μ L, ecopipette	12.00 / 2.00	6.00 / 1.00	●
C02-1	Capp pipette, variable vol. 0.2-2 μ L, ecopipette (length as C10-1)	12.00 / 2.00	6.00 / 1.00	●
C10-1	Capp pipette, variable vol. 0.5-10 μ L, ecopipette	2.50 / 1.00	1.80 / 0.50	○
C20-1	Capp pipette, variable vol. 2-20 μ L, ecopipette	5.00 / 1.00	1.50 / 0.30	●
C50-1	Capp pipette, variable vol. 5-50 μ L, ecopipette	2.00 / 1.00	1.20 / 0.30	●
C100-1	Capp pipette, variable vol. 10-100 μ L, ecopipette	2.50 / 0.70	1.00 / 0.20	●
C200-1	Capp pipette, variable vol. 20-200 μ L, ecopipette	2.50 / 0.60	0.70 / 0.20	●
C1000-1	Capp pipette, variable vol. 100-1000 μ L, ecopipette	2.00 / 0.80	0.60 / 0.20	●
C5000-1	Capp pipette, variable vol. 1-5mL, ecopipette	1.20 / 0.60	0.30 / 0.15	●

ecopipette fixed volume

Cat. No.	Description	Inaccuracy %	Imprecision %
C02-1S-001	Capp pipette, Fixed Volume 1 μ L , ecopipette	2.50	1.80
C10-1-005	Capp pipette, Fixed Volume 5 μ L, ecopipette	1.20	0.80
C10-1-010	Capp pipette, Fixed Volume 10 μ L, ecopipette	1.00	0.50
C50-1-020	Capp pipette, Fixed Volume 20 μ L, ecopipette	1.00	0.30
C50-1-025	Capp pipette, Fixed Volume 25 μ L, ecopipette	0.80	0.35
C100-1-030	Capp pipette, Fixed Volume 30 μ L, ecopipette	1.16	0.40
C100-1-050	Capp pipette, Fixed Volume 50 μ L, ecopipette	1.00	0.30
C200-1-100	Capp pipette, Fixed Volume 100 μ L , ecopipette	0.70	0.20
C200-1-200	Capp pipette, Fixed Volume 200 μ L, ecopipette	0.60	0.20
C1000-1-500	Capp pipette, Fixed Volume 500 μ L, ecopipette	0.80	0.15
C1000-1-1000	Capp pipette, Fixed Volume 1000 μ L, ecopipette	0.80	0.20
C5000-1A-5000	Capp pipette, Fixed Volume 5000 μ L, comfort	0.60	0.15
CXXX-1-YYY	Other fixed vol. pipettes (xxx=max pipette volume; yyy=desired volume)		