

LAB MONITORING & ALARM SYSTEM



» Experts In Monitoring And Alarm Systems



» Text Message Alarms to Mobile Phone



» Powerful Wireless Signals

Monitoring System Britannia

PROVEN SOLUTIONS FOR ALL IVF EQUIPMENT MONITORING

» Completely Failsafe

MONITORING STAND-ALONE UNIT

A stand-alone sensing unit with easy user recalibration function, remote flashing xenon lamp and sounder includes local display to show current reading of gas concentration.



Carbsafe CO2 Monitor



Oxysafe O2 Depletion Monitors

Transmitters

Single and multi-channel transmitters capable of measuring and transmitting Temperature

- Humidity
- Voltage/Current
- Visible light/UV
- Event
- Differential pressure
- Gas ; CO2/O2



DARCA H1 Software

Suite enables users to configure, download and analyses the data transmitted and collected

Sensors



Temperature, Humidity, CO2, Differential pressure and many other environmental sensors, ideally suited to the most demanding of applications



Main Panel Base Station



ALARM TEXT MESSAGE



MOXA Communicate with the system over the LAN from multiple PC's

Labsafe Key Points



250 Channels

Ability to monitor up to 250 channels, input or sensor



Transmitted data is encrypted



2 km Transmissions

System can transmit up to 2 km over open ground



SMS Alarm

Can text up to 10 phone numbers



Battery backup