

AFIAS SERIES

The AFIAS Series represents a fluorescence immunoassay-based POC analyzer, featuring an integrated Sample-to-Answer system that can conduct 1, 3, and 6 parallel tests based on the model. This compact POC system delivers prompt and reliable results in real-time. The AFIAS All-in-one cartridge further simplifies the process, enhancing user convenience and facilitating seamless operation.

Automated Fluorescence Immuno Assay System



AFIAS-1

FEATURES

- Simple procedure with All-in-One cartridges
- Easy to use without pipetting with proprietary C-tip
- Quick test results with reliability
- Small, compact, & light

AFIAS-3

FEATURES

- AFIAS system for 3 parallel tests
- Extra tip station – faster TAT
- Compact system for multiple testing
- User friendly UI



AFIAS-6

FEATURES

- Sample-to-Answer POCT system
- Easy to use without pipetting with proprietary C-tip
- Up to 6 different parallel tests
- Reliable test results

