


# QUALITY CONTROL

## NILOCHECKER 500



### One Reference instrument for many applications

With the NiloChecker you can control and calibrate incubators, laminar flow workstations, warming system etc. No need for multiple instruments from various manufacturers. An increasing range of measurement probes ensures your investment in years to come.

#### FEATURES

- The all-in-one solution
- Measure CO2 and O2 and temperature inside any incubator
- Save precious time on quality and process controls
- High battery capacity
- No maintenance or calibration

SPECIFICATIONS	
Measure range and accuracy	Depends on connected probe
Display and Touch	7", 800 x 480 pixel Adjustable brightness Capacitive touch
Compliance	ISO/EN 61010:2010 / IEC 60950-1:2005 / AMD2:2013 / CE
Operating conditions	0-50 °C, 5-95% RH, IP 50
Power supply	12 V 30 W, Wall plug adaptor, 100-240V
Connectors	5 Pcs Stainless steel with lock and quick release
Battery time between charging	8 Hours Level displayed in %
Battery Charging Time	< 3 hours

#### DIMENSIONS

