


WORKSTATIONS

EnviSafe *uv*

designed for optimal user and sample protection.

It is prepared for a stereo microscope for observation under a cabinet. It comes with built-in heated stage under the microscope to provide optimal working condition.




FEATURES

- Multiple protection design increases security level significantly
- Innovative surrounded air intake grill design
- 1,000 Lux fluorescent lamp at the work surface
- Germicidal UV light system
- Ergonomic 10 degrees angle window design
- Equipped with exhaust HEPA filter 99.99% efficiency @ 0.3 micron (Class H14)
- 7" High resolution touch screen display
- Mounting assembly for most stereo microscopes integrated with strengthened heated glass
- Table top is made of stainless steel integrated with 2-sides heating areas to give an optimal working condition



TOUCH SCREEN DISPLAY



**TOKAI HIT THERMO PLATE
TPIE-110RX (OPTIONAL)**

