

INCUBATORS

Synvivo™


ensures the optimal culture environment for the best choice among right embryos.

FEATURES

- Synvivo is like a little apartment complex composed of eight individual independently working mini-incubators
- Its eight chambers have individual on/off control for their gas supply
- Each & individual chamber is equipped with heating elements on top as well as bottom of the chambers
- 7" touch screen offers exact display configuration for user's convenience
- A VOC Filter administration program is integrated to automatically alert filter replacement
- Compatible with most of the culture dishes
- Easy administration of log data is possible through a USB port
- Quality certified components with proven performance allow Synvivo to produce optimal culture environment



Benchtop Multi-chamber Incubator



Every dish is seated securely and comfortably

- 16 x Nunc 4well dishes
- 16 x Falcon / Nunc 60mm dishes
- 16 x 45mm multi well dishes
- 48 x Nunc 35mm dishes

